

## WATER BATH

Product Code : EEWB10

HSN Code : 841990

Thermostatically Water Bath is a double walled chamber made of SS and outer made of MS with powder coating. The gap between the inner and outer chamber fines with special grade glass-wool insulation to avoid the heat loss. Concentric rigs can be removed to accommodate containers of various sizes.



## PRODUCT SPECIFICATION

Temp. Range	5 degree C above ambient to 99 degree C
Brand	Elena
Accuracy	1 degree C
Power Supply	230 Volts A.C 50 Hz
Phase	Single Phase
Body Material	Stainless Steel
Surface	Mild Steel powder coated

## AVAILABLE SIZES:

Size	Rating Watts	Holes
35 x 25 x 10 cm	1000	6 Hole of 3"Dia
40 x 30 x 10 cm	1500	12 Hole of 3"Dia